## A peer-reviewed version of this preprint was published in PeerJ on 6 May 2014.

<u>View the peer-reviewed version</u> (peerj.com/articles/377), which is the preferred citable publication unless you specifically need to cite this preprint.

Craft WF, Conway JA, Dark MJ. 2014. Comparison of histomorphology and DNA preservation produced by fixatives in the veterinary diagnostic laboratory setting. PeerJ 2:e377 <u>https://doi.org/10.7717/peerj.377</u>

# Comparison of histomorphology and DNA preservation produced by fixatives in the veterinary diagnostic laboratory setting

**Abstract:** Histopathology is the most useful tool for diagnosis of a number of diseases, especially cancer. To be effective, histopathology requires that tissues be fixed prior to processing. Formalin is currently the most common histologic fixative, offering many advantages: it is cheap, readily available, and pathologists are routinely trained to examine tissues fixed in formalin. However, formalin fixation substantially degrades tissue DNA, hindering subsequent use in diagnostics and research. We therefore evaluated three alternative fixatives, TissueTek® Xpress® Molecular Fixative, modified methacarn, and PAXgene®, all of which have been proposed as formalin alternatives, to determine their suitability for routine use in a veterinary diagnostic laboratory.

This was accomplished by examining the histomorphology of sections produced from fixed tissues as well as the ability to amplify fragments from extracted DNA. Tissues were sampled from two dogs and four cats, fixed for 24-48 hours, and processed routinely. While all fixatives produced acceptable histomorphology, formalin had significantly better morphologic characteristics than the other three fixatives. Alternative fixatives generally had better DNA amplification than formalin, although results varied somewhat depending on the tissue examined. While no fixative is yet ready to replace formalin, the alternative fixatives examined may be useful as adjuncts to formalin in diagnostic practices.

1	William F. Craft <sup>1</sup> , Julia A. Conway <sup>1</sup> , and Michael J. Dark <sup>1,2</sup>		
2			
3	1 – Department of Infectious Diseases and Pathology, College of Veterinary		
4	Medicine, University of Florida, Gainesville, FL		
5	2 – Emerging Pathogens Institute, University of Florida, Gainesville, FL		
6			
7	Corresponding Author:		
8	Michael J. Dark		
9	Department of Infectious Diseases and Pathology		
10	College of Veterinary Medicine		
11	P.O. Box 110880		
12	Gainesville, FL 32611-0880		
13	Phone: (352) 294-4138		
14	E-mail: <u>darkmich@ufl.edu</u>		
15			

#### 16 Introduction:

17 Histopathology is the most useful tool for diagnosis of a number of diseases, 18 especially cancer. To be effective, histopathology requires that tissues be fixed prior to processing. The ultimate tissue fixative for histopathology would create tissue 19 20 histomorphology identical to formalin, pose no hazard to human health, preserve nucleic material for an extended period of time preferably at room temperatures, and be cost 21 effective. Formalin has been the most common histologic fixative for over 100 years. It 22 offers many advantages: it is cheap, readily available, and pathologists are routinely 23 24 trained to examine formalin-fixed tissues (Gugic et al. 2007; Srinivasan et al. 2002).

However, formalin has several disadvantages as well. It can cause respiratory
irritation and is classified as a carcinogen (Bolt et al. 2010; Bosetti et al. 2008; Buesa
2008; Gugic et al. 2007). It can also degrade nucleic acids and proteins, which make
formalin-fixed tissues less usable for downstream molecular diagnostics (Buesa 2008;
Gugic et al. 2007). Formalin is an aldehyde-based fixative that works by cross-linking
proteins, which irreversibly degrades proteins and nucleic acids (Srinivasan et al. 2002).

The most common method for preserving tissues for molecular diagnostics is freezing at -80°C, although this does not allow for histopathologic examination of tissues. This method itself has a number of issues, including requiring special equipment, difficulty in shipping samples to laboratories, and requiring duplicate samples to be taken for histopathologic examination.

This is impractical for many private veterinary practices, as these generally lack 36 37 the facilities to freeze samples at -80°C and many samples are too small to duplicate samples. Therefore, alternative fixatives have been proposed to allow for both 38 39 histopathologic examination and molecular diagnostics {Cox, 2006 #2;Gugic, 2007 #3; Kap, 2011 #6; Vincek, 2003 #10}. These have been shown to preserve nucleic acids 40 with results similar to those obtained with fresh or frozen tissues, while still preserving 41 histomorphology. Some of the more successful alternatives include Tissue-Tek® 42 43 Xpress® Molecular Fixative (Gugic et al. 2007), PAXgene® (Kap et al. 2011), and modified methacarn solution (Cox et al. 2006). All of these are alcohol-based and non-44 cross-linking. Evaluation of histomorphology preservation has varied among the studies 45 evaluating these fixatives and generally involve research settings using techniques that 46 are not practical in most clinical situations. In addition, some of these fixatives are 47 currently cost prohibitive in the veterinary clinical setting. The majority of the studies 48 evaluating these fixatives have evaluated single organs from humans or rodents {Cox, 49 2006 #2} or multiple organs from humans {Kap, 2011 #6; Vincek, 2003 #10}. One study 50 evaluated Tissue-Tek® Xpress® Molecular Fixative and formalin comparing 51 histomorphology and RNA quality from a variety of animal tissues (small animals, 52 rodents, lagamorphs, birds, insects, and lizards) both at room temperature and high 53 ambient temperatures simulating field collection of samples {Gugic, 2007 #3}. They 54 concluded that Tissue-Tek® Xpress® Molecular Fixative protected RNA and provided 55 56 acceptable histomorphology that would not hinder histologic diagnosis in the species studied. Some studies evaluating multiple animal species have included limited numbers 57 of fixatives for comparison {Gugic, 2007 #3; Vincek, 2003 #10}. 58

The main limitation of all of these previous studies is that they have evaluated 59 fixatives in a research setting. There has not been a systematic evaluation of these to 60 determine their utility in the veterinary diagnostic setting. Alternative fixatives would 61 have a number of benefits for veterinary diagnostic laboratories, including (depending on 62 the nature of the fixative) decreasing hazardous waste disposal costs, decreasing health 63 risks to laboratory workers, and enhancing the power of retrospective studies Therefore, 64 we conducted this study to determine how alternative fixatives would function in a 65 standard diagnostic laboratory setting by evaluating histomorphology of a variety of 66 tissues from dogs and cats, as well performing a quantitative evaluation of recoverable 67 DNA from tissues. 68

69

#### 70 Materials and Methods:

All study protocols were approved by the University of Florida Institutional 71 Animal Use and Care Committee (approval #201105654), and all animals were 72 euthanized for reasons unrelated to this project. Necropsies were performed on four cats 73 74 and two dogs within four hours of euthanasia. Replicate 10 x 10 x 5 mm samples from the liver, brain, lung, lymph node, kidney and spleen were collected. One sample of each 75 76 tissue was frozen at -80°C. The remaining samples were placed into 10% neutral buffered 77 formalin, Tissue-Tek® Xpress® Molecular Fixative (TT-XMF), modified methacarn, and PAXgene<sup>®</sup>, with a minimum of 1:10 tissue to fixative volume. Samples were allowed to 78 fix for 24-48 hours at room temperature with the exception of tissues in PAXgene®, 79 which were fixed and preserved according to the manufacturer's protocol. 80

All tissues were processed using a Tissue-Tek processor with a standard overnight protocol (excluding formalin steps) followed by paraffin embedding and hematoxylin and eosin staining. The 10% neutral buffered formalin, Tissue-Tek® Xpress® Molecular Fixative, and PAXgene® fixatives were purchased commercially (ThermoFisher Scientific, Waltham, MA). Modified methacarn was prepared as previously described, using 8 parts methanol and 1 part glacial acetic acid (Cox et al. 2006).

Histomorphology was evaluated by two blinded board-certified veterinary 87 anatomic pathologists (MJD, JAC) and one blinded anatomic pathology resident (WFC). 88 Histomorphology of nuclear, cytoplasmic, and cellular membrane detail were evaluated 89 on a 1-4 scale (table 1). Sample scores were averaged between all three evaluators. For 90 one cat, the formalin-fixed lymph node sample was lost from the block; therefore, 91 formalin fixation histomorphometry scores for lymph node are based on the remaining 92 93 five samples. Both the individual components of the histomorphometry score as well as the total score were evaluated using a Kruskal-Wallis test (Lowry 2012) to determine if 94 there was a difference between any of the four groups. If a significant difference was 95 found (p < 0.05), the Mann-Whitney test was used to compare each group to each other 96

group, to determine significant differences between each individual fixative. Fixatives
were considered significantly different if the one-tailed Mann Whitney p value was less
than 0.05. The minimum, 25<sup>th</sup> quartile, median, 75<sup>th</sup> quartile, and maximum were
calculated for each tissue as well as for all tissues combined using Microsoft Excel
(v14.3.9, Microsoft Corp., Seattle, WA). Graphs were generated using GNUplot (v.4.6,
patchlevel 3).

Tissue scrolls were obtained from the paraffin blocks one week after processing 103 and DNA was extracted using the QIAamp DNA FFPE Tissue kit (Qiagen Inc., Valencia, 104 CA). As a control, DNA was extracted from tissues frozen at -80°C using the OIAamp 105 DNA Mini Kit (Oiagen Inc.). Primers were designed by aligning the sequences of the 106 retinol-binding protein 3, interstitial gene (IRBP) from dog, mouse, rat, and human, and 107 selecting regions that were relatively conserved, to generate 100, 200, 300, 500 and 750 108 base pair long amplicons (table 2). Extracted DNA was amplified via PCR on an Applied 109 Biosystems Veriti Thermal Cycler with the following conditions: 96°C for 3 minutes, 110 followed by 35 cycles of 96°C for 1 minute, 60°C for 1 minute, then 72°C for 1 minute. 111 This was followed by 7 minutes at 72°C, with a final hold at 4°C until the next morning. 112 Samples were examined on a 1.25% agarose gel via electrophoresis. 113

The presence or absence of bands for all sizes was noted. The Kruskal-Wallis test was used to determine if there was a significant difference between the maximum band size for any of the fixatives. If significant (p<0.05), the Mann-Whitney test was used to compare each fixative against each other to determine which had significant differences (one-tailed p<0.05). The minimum,  $25^{th}$  quartile, median,  $75^{th}$  quartile, and maximum calculated for each fixative and for each tissue using Microsoft Excel. Graphs were generated using GNUplot (v.4.6, patchlevel 3).

#### 121 **Results:**

#### 122 Histomorphology

While the majority of the alternative fixatives produced adequate histomorphology in the tissues examined, formalin fixed tissues consistently resulted in superior histomorphology. There was no statistically significant difference between mean histomorphology scores comparing dog and cat tissues, and these were combined for subsequent analysis. Mean total, nuclear, cytoplasmic, and cellular membrane scores (figs. 1A-D) for formalin fixed tissues were higher than for all other fixatives (p<0.0001), although there is substantial variation with all fixatives (figs. 2A-D).

While a number of minor artifacts were noted, the primary difference notedbetween formalin and the other fixatives was in erythrocytes. This is likely reflected in

the significantly higher scores for formalin vs. other fixatives in the spleen (p=0.0026),an organ made up in large part by erythrocytes.

#### 134 DNA Preservation

Formalin has significantly shorter total maximum DNA band sizes than TT-XMF (p=0.0158), modified methacarn solution (p < 0.0001), and PAXgene (p=0.0004) (fig. 3). In particular, the bands obtained from lymph nodes were significantly smaller with formalin than with TT-XMF (p=0.0179) and modified methacarn (p=0.004) (fig. 4). Overall, modified methacarn solution performed as well or better than the other fixatives for all tissues, with the best score in brain (median amplicon length of 750bp).

### 141 **Discussion:**

While alternative fixatives have been found to work well in research settings (Cox 142 et al. 2006; Kap et al. 2011; Vincek et al. 2003), these are not ready to replace formalin 143 for routine tissue processing in the veterinary laboratory. All of the fixatives require 144 tissues be prevented from contacting formalin to benefit from their nucleic acid 145 preserving qualities, which would require laboratories to either maintain separate tissue 146 processors or bar submission of formalin-fixed tissues. Neither of these is practical in 147 veterinary practice. Several fixatives produce excellent histomorphology with alternative 148 processing techniques; this is also impracticable in most veterinary diagnostic 149 150 laboratories, as it would require separate processing runs.

However, while no fixative is ideal from the standpoint of replacing formalin, all 151 fixatives produced interpretable slides. Therefore, using alternative fixatives may be 152 153 useful in specific circumstances where subsequent DNA isolation may be required. For example, tissue samples from neoplasms may be saved separately to generate a tissue 154 bank for subsequent research projects. The specific alternative chosen should be based on 155 the tissue selected, as well as predicted needs for DNA amplification and preservation of 156 histomorphology. For example, while TT-XMF had better histomorphology scores in the 157 kidney than either modified methacarn or PAXgene, it had a lower median DNA 158 amplicon size. 159

One characteristic observed with alternative fixatives was that bloody or 160 161 congested tissues often had unfixed areas, which could result in missing lesions and inaccurate diagnoses. This has not been found in previous studies (Cox et al. 2006), and 162 may be due to a number of factors. First, the size of sample taken will greatly influence 163 fixation. For most veterinary diagnostic laboratories, 1 cm thick samples are considered 164 standard for histopathologic examination. In many previous studies, samples taken for 165 fixation were substantially thinner; for example, the study by Cox et al. used 15 mm x 8 166 mm x 3 mm samples. Other possibilities include differences in processing; microwave 167

fixation (Cox et al. 2006)} or rapid tissue processing (Vincek et al. 2003) techniques have
been used. Tuning the processing technique for the fixative selected would likely
improve fixation and the ultimate histomorphology.

Finally, our evaluation of macromolecule preservation of was limited to DNA. 171 Additional analysis would be required to determine whether these fixatives preserve 172 RNA equally well. Other variables require investigation to determine the best fixative for 173 a particular application. These include the effects of fixation time on nucleic acid quality, 174 as many samples will sit longer than 24 hours before processing, as well as the effect of 175 storage time after tissue processing but before sectioning for nucleic acid isolation, since 176 177 many blocks will be stored for a period of time between the evaluation of histopathology and nucleic acid isolation. The latter is especially important if laboratories set up tissue 178 banks, as samples would be expected to be stored for prolonged periods. 179

#### **180 Conclusions:**

While no fixative is ideal to replace formalin, alternative fixatives have generally acceptable histomorphologic characteristics in most tissues and are valuable adjuncts to standard formalin fixation. Projects proposing to use an alternative fixative for a research project should first evaluate the project requirements and ideally test the fixative with samples of the target organ to determine the best fixative, processing characteristics, and histomorphology compromises before actual sample collection begins.

#### 187 Acknowledgements:

The authors would like to acknowledge Patrick Knisley for his assistance with samples,
the University of Florida College of Veterinary Medicine Histopathology Laboratory for
help with tissue processing, and Antoinette McIntosh for help with DNA processing.

191 **References:** 

192

193	Bolt HM, Degen GH, and Hengstler JG. 2010. The carcinogenicity debate on
194	formaldehyde: How to derive safe exposure limits? Archives of toxicology
195	84:421-422.
196	Bosetti C, McLaughlin JK, Tarone RE, Pira E, and La Vecchia C. 2008. Formaldehyde
197	and cancer risk: a quantitative review of cohort studies through 2006. Annals of
198	oncology : official journal of the European Society for Medical Oncology / ESMO
199	19:29-43.
200	Buesa RJ. 2008. Histology without formalin? Annals of diagnostic pathology 12:387-396.
201	Cox ML, Schray CL, Luster CN, Stewart ZS, Korytko PJ, KN MK, Paulauskis JD, and
202	Dunstan RW. 2006. Assessment of fixatives, fixation, and tissue processing on
203	morphology and RNA integrity. Experimental and molecular pathology 80:183-
204	191

205	Gugic D, Nassiri M, Nadji M, Morales A, and Vincek V. 2007. Novel tissue preservative
206	and tissue fixative for comparative pathology and animal research. Journal of
207	Experimental Animal Science 43:271-281.
208	Kap M, Smedts F, Oosterhuis W, Winther R, Christensen N, Reischauer B, Viertler C,
209	Groelz D, Becker KF, Zatloukal K et al 2011. Histological assessment of
210	PAXgene tissue fixation and stabilization reagents. <i>PloS one</i> 6:e27704.
211	Lowry R. 2012. VassarStats: Website for Statistical Computation (accessed 3/1/2012
212	2012).
213	Srinivasan M, Sedmak D, and Jewell S. 2002. Effect of fixatives and tissue processing on
214	the content and integrity of nucleic acids. The American journal of pathology
215	161:1961-1971.
216	Vincek V, Nassiri M, Nadji M, and Morales AR. 2003. A tissue fixative that protects
217	macromolecules (DNA, RNA, and protein) and histomorphology in clinical
218	samples. Laboratory investigation; a journal of technical methods and pathology
219	83:1427-1435.
220	
221	Figure Legends:
222	Figure 1. Histomorphology scores for all animals and tissues combined. The median is
223	represented by a red diamond, the box represents the 25 <sup>th</sup> and 75 <sup>th</sup> quartiles, and the
224	whiskers represent 1.5 x interquartile range
	ministers represent i.e a interquitire funde.

Figure 2. Histomorphology scores for individual tissues. The median is represented by a red diamond, the box represents the  $25^{th}$  and  $75^{th}$  quartiles, and the whiskers represent 1.5 x interquartile range.

Figure 3. DNA amplicon size ranges for all animals and tissues combined. The median is represented by a red diamond, the box represents the 25<sup>th</sup> and 75<sup>th</sup> quartiles, and the whiskers represent 1.5 x interquartile range.

Figure 4. DNA amplicon sizes for different tissue samples. The median is represented by

a red diamond, the box represents the  $25^{\text{th}}$  and  $75^{\text{th}}$  quartiles, and the whiskers represent 1.5 x interquartile range.

Figure 5. Representative fixative histomorphology. Samples are from the liver of a single
cat. A – formalin, B – TT-XMF®, C – modified methacarn, D – PAXgene®.

236 Table Legends:

237 Table 1. Histomorphology scoring chacteristics.

Table 2. DNA primers used in this study.

239









244





250

**PeerJ** PrePrints



252 Fig 5.

253

**PeerJ** PrePrints

# 254

Characteristic	Score	Criteria
Nuclear	4	Sharp nuclear membrane; chromatin pattern and nucleolus, when present, are distinct
	3	Slight degradation in chromatin pattern
		Nucleolus, when present, less distinct but discernable, sharp nuclear membrane
	2	Less distinct nuclear membrane, fuzzy chromatin pattern
		Nucleolus, when present, is difficult to discern
	1	Fuzzy nuclear membrane
		Chromatin pattern difficult to determine
		Nucleoli indetectable
	0	Nucleus not able to be differentiated from cytoplasm
Cytoplasm	4	Normal cellular morphology easily determined
	3	Intracytoplasmic details fuzzy
	2	Only rare evidence of normal intracellular structures
	1	Increased cytoplasmic pallor, increased cytoplasmic eosinophilia
	0	Cytoplasm homogenously pale eosinophilic with no evidence of organelles
	4	Cells have distinct intracellular borders
Cell Membranes		Normal substructures, if present, are easily distinguished
	3	Loss of substructures in some cells
		Slight loss of intracellular details
	2	Loss of substructures in most cells
		Obvious fuzzing of many cellular borders
	1	No substructures detected
		Significant fuzzing of most cellular borders
TD 11	0	Cells unable to be distinguished from adjacent cells
255 Table	1	
256		

Primer Name	Sequence
IRBP_F	CCT KGT RCT GGA NAT GGC
IRBP_R1_100bp	CTC TTG ATG GCC TGC TC
IRBP_R2_200bp	GGC TCA TAG GAG ATG ACC AG
IRBP_R3_300bp	CAG GTA GCC CAC RTT NCC CTC
IRBP_R4_400bp	CGG AGR TCY AGC ACC AAG G
IRBP_R5_500bp	GAT CTC WGT GGT NGT GTT GG
IRBP_R6_750bp	CTC AGC TTC TGG AGG TCC

258 Table 2